

Portable Oxygen Concentrators • FAA Approved • Battery Duration



Weight: 8.5 lb  
Dims: 12" x 6" x 8.5"

**Phillips Respironics EverGo™**

1 Battery / with 2nd Battery  
Setting 1 = 6 Hr / 12 Hr  
Setting 2 = 4 Hr / 8 Hr  
Setting 3 = 3 Hr / 6 Hr  
Setting 4 = 2 Hr / 4 Hr  
Setting 5 = 2 Hr / 4 Hr  
Setting 6 = 2 Hr / 4 Hr



Weight: 9.7 lb  
Dims: 11.62" x 6" x 12.39"

**Inogen One™**

Setting 1 = 3.0 Hr  
Setting 2 = 3.0 Hr  
Setting 3 = 2.5 Hr  
Setting 4 = 2.5 Hr  
Setting 5 = 2.0 Hr



Weight: 4.9 lb  
Dims: 9.5" x 7.5" x 3.125"

**Inova Labs International Biophysics™ LifeChoice**

Internal / with External Battery  
Setting 1 = 2 Hr / 5 Hr  
Setting 2 = 2 Hr / 5 Hr  
Setting 3 = 2 Hr / 5 Hr



Weight: 6 lb  
Dims: 10" x 7" x 4"

**Invacare XPO2 - XPO100™**

Internal / with External Battery  
Setting 1 = 3.5 Hr / 7 Hr  
Setting 2 = 2.5 Hr / 5 Hr  
Setting 3 = 2.0 Hr / 5 Hr  
Setting 4 = 1.5 Hr / 3 Hr  
Setting 5 = 1.0 Hr / 2 Hr



Weight: 19 lb  
Dims: 15" x 11" x 8"

**DeVilbiss Healthcare iGo™**

Pulse / Continuous  
Setting 1 = 5.4 Hr / 4.0 Hr  
Setting 2 = 4.7 Hr / 2.4 Hr  
Setting 3 = 4.0 Hr / 1.6 Hr  
Setting 4 = 3.5 Hr  
Setting 5 = 3.2 Hr  
Setting 6 = 3.0 Hr



Weight: 14.85 lb  
Dims: 8" x 8" x 12"

**Oxygo2 LLC Independence™**

Pulse / Continuous  
Setting 1 = 2.8 Hr / 2.4 Hr  
Extended Battery = 5.6 Hr / 4.8 Hr  
Setting 2 = 2.8 Hr / 1.4 Hr  
Extended Battery = 5.6 Hr / 2.8 Hr  
Setting 3 = 1.4 Hr / 53 min  
Extended Battery = 2.8 Hr / 1.8 Hr  
Setting 4 = 1.4 Hr  
Extended Battery = 2.8 Hr  
Setting 5 = 53 min  
Extended Battery = 1.8 Hr  
Setting 6 = 53 min  
Extended Battery = 1.8 Hr



Weight: 4.4 lb  
Dims: 8.6" x 6.1" x 3.6"

**AirSep FreeStyle™**

Internal Battery / with Battery Belt  
Setting 1 = 3.5 Hr / 10 Hr  
Setting 2 = 2.5 Hr / 6 Hr  
Setting 3 = 2.0 Hr / 5 Hr



Weight: 9.75 lb  
Dims: 5.5" x 7.25" x 16.31"

**AirSep LifeStyle™**

Setting 1 = 3.25 Hr\*  
Setting 2 = 3.25 Hr\*  
Setting 3 = 3.25 Hr\*  
Setting 4 = 3.25 Hr\*  
Setting 5 = 3.25 Hr\*  
\*with Power Pack (4 batteries)  
Each Battery lasts 50 min



Weight: 9.8 lb  
Dims: 11.6" x 7.4" x 4.6"

**Delphi Central Air, Evo™**

Setting 1 = 3 Hr 25 min  
Setting 2 = 2 Hr 45 min  
Setting 3 = 2 Hr 15 min  
Setting 4 = 1 Hr 40 min  
Setting 5 = 1 Hr 15 min

**Editor's Note: At press time, this Delphi unit is currently unavailable.**



Weight: 7 lb  
Dims: 11.61" x 5.98" x 10.75"

**Inogen One G2™**

Standard / Extended  
Setting 1 = 4.0 Hr / 8.0 Hr  
Setting 2 = 3.5 Hr / 7.0 Hr  
Setting 3 = 2.7 Hr / 5.4 Hr  
Setting 4 = 2.2 Hr / 4.4 Hr  
Setting 5 = 2.0 Hr / 4.0 Hr



**Invacare Solo™**

Pulse / Continuous  
Setting .5 = N/A / 4.5 Hr  
Setting 1 = 4.5 Hr / 3.5 Hr  
Setting 2 = 3.5 Hr / 2.5 Hr  
Setting 3 = 3.0 Hr / 1.5 Hr  
Setting 4 = 2.5 Hr  
Setting 5 = 2.5 Hr



Weight: 17.9 lb  
Dims: 19.3" x 12.3" x 7.1"

**SeQual Eclipse®**

Pulse / Continuous  
Setting .5 = NA / 4.4 Hr  
Setting 1 = 4.4 Hr / 3.7 Hr  
Setting 2 = 3.6 Hr / 2.4 Hr  
Setting 3 = 3.0 Hr / 1.3 Hr  
Setting 4 = 2.6 Hr  
Setting 5 = 2.3 Hr  
Setting 6 = 2.1 Hr

Battery Duration Chart Estimated Times Based on NEW Batteries  
Hr = hours / min = minutes / lb = pounds / Dims = dimensions



POC information courtesy of Trevor Scribner  
Advanced Aeromedical, Inc.  
www.aeromed.com